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Deutschland

Contents

| | | |
|----|-------------------------------|----|
| 1. | Important instructions..... | 6 |
| 2. | Included in delivery | 7 |
| 3. | Product features..... | 8 |
| 4. | Operating..... | 8 |
| 5. | Cleaning and maintenance..... | 10 |
| 6. | Disposal | 10 |
| 7. | Warranty..... | 11 |

1. Important instructions

- Non-observance of the following notes may result in personal injury or material damage.
- This heated belt is only designed to warm up the human body.
- Do not use the heat belt on children or persons who are disabled, asleep or insensitive to heat (persons who cannot react on overheating).
- Children shall not play with the appliance.
- Do not pull, twist or make sharp bends in the cables.
- If the heating belt is being used over several hours, we recommend that you set the lowest temperature setting on the battery to avoid overheating the heated body part, which may lead to burns to the skin.

2. Included in delivery

- 1 x heated belt
- 1 x USB charger
- 1x battery bag
- 1 x Extension USB cable
- 1 x Extension for belt
- 1 x Instruction Manual

3. Product features

The Thermrup **Far infrared Wireless Heated Belt** is designed for indoor and outdoor use provides deep heat with far infrared for back and abdomen with good mobility. It is easy to wear, comfortable and safe. It has 3 adjustable temperature settings. The special rubber heating element has large heating area and even heat distribution. The belt can be washed in washing machine.

4. Operation

4.1 Operating with power bank

- Place the heated belt around your back and abdomen and secure the belt by fastening the Velcro. Depending on the position you put the belt on, it warms different areas of your body
- Plug your power bank (not included) into the battery bag.
- Plug the USB end of the cable into the USB port of the power bank. Turn on the power bank.
- Puts the power bank in the battery bag. Place the battery bag on the heated belt.
- Press the power button of the controller for approx. 2 seconds to switch on the controller. Then press the power button briefly (if necessary several times) to set the desired temperature level (1-2-3, red-blue-green).
The 1st stage is the warmest stage, the LED indicator lights up red.
The 2nd stage is the middle stage, the LED indicator lights up blue.
The 3rd stage is the low level, the LED indicator is green.
- Press the power button for approx. 2 seconds to switch off the controller.

4.2 mains operation

- Connect the extension cable to the connection cable of the belt.
- Insert the USB plug of the extension cable into the USB charger.
- Press the power button of the controller for approx. 2 seconds to switch on the controller. Then press the power button briefly (if necessary several times) to set the desired temperature level (1-2-3, red-blue-green).
- Press the power button for approx. 2 seconds to switch off the controller.
- After use, remove the extension cord from the charger.

*** Please use only the supplied charger. ***

5. Cleaning and maintenance

- Before cleaning, disconnect the heat belts from the power bank or power supply. Take the heating pad and the cable out of the belt.
- Never use aggressive cleaning agents or strong brushes.
- The heat belt can be washed in the washing machine up to 45 ° C in a washing bag or supplied washing bag.
- The thermal belt must not be dried or ironed in a dryer.

6. Disposal

This product must not be disposed together with the domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

7. Warranty

We provide a 2-year warranty from the date of purchase for material and manufacturing defects of the product.

The warranty does not apply:

- In the case of damage caused by improper use.
- For wear parts.
- Customer's own fault.
- Once the unit has been opened by an unauthorized workshop.

For the assertion of a warranty within the guarantee period, the customer must prove the purchase.

In case of complaints please contact our service under the following contact:

E-Mail: service@zewox.de

Zewox GmbH

Oststr.159

47057 Duisburg

Germany